



Efficient Workflow Precision Results

EasyFrost slides and EasyFrost adhesion slides from EpreDia
Corrosion resistant glass microscope slides with ground edges

For standard laboratory applications:

- Histology
- Hematology
- Microbiology
- Urinalysis





EasyFrost Adhesion Microscope Slides

EasyFrost adhesion microscope slides feature a positive charge coating, allowing firm attraction to tissue samples. Manufactured to high quality standards in our Shanghai manufacturing site, these slides are made of highly transparent white glass.

EasyFrost Microscope Slides

EasyFrost microscope slides are manufactured to high quality standards in our Shanghai manufacturing site, and made of corrosion resistant green glass. These slides are suitable for routine laboratory applications.



- All slides are pre-cleaned with ground edges for smoother feel
- Corrosion resistant glass packed 50 slides per box
- Each slide features a white frosted end for easy-labeling
- Each pack is individually wrapped to protect against moisture and humidity
- EasyFrost adhesion slides feature an adhesive coating for clinical applications
- Option includes clipped corners to reduce danger of cuts

EasyFrost Adhesion Microscope Slides

- White frosted end, 25 x 75 x 1 mm, ground edge 90°
- 50 slides/pack, 40 packs/case 32550T
- Manufactured from highly transparent white glass

- White frosted end, 25 x 75 x 1 mm, ground edge 90°, clipped corners
- 50 slides/pack, 40 packs/case 8519W-ADC-4
- Manufactured from highly transparent white glass

EasyFrost Microscope Slides

- White frosted end, 25 x 75 x 1 mm, ground edge 90°, clipped corners
- 50 slides/pack, 40 packs/case 8106W-4
- Manufactured from corrosion resistant green glass

EpreDia Laboratory Products Manufacturing (Shanghai) Co., Ltd
T6-7, No. 5001 Huadong Road
Jianqiao, Shanghai, China
Tel: +86 021 58307022
Fax: +86 021 58307022
Email: ep-sales.cn@epredia.com

Find out more at epredia.com

© 2020 EpreDia Group. All rights reserved. All trademarks are the property of EpreDia Group and its subsidiaries unless otherwise specified. These products are not available in all regions. Please reach out to your EpreDia sales representative for more information. **M41058 R112020**

 Enhancing precision
cancer diagnostics

POWERED BY SHANDON MICRON Menzel-Gläser